# Cuvette washer


## Operating manual for cuvette washers 2 990 111 (211)

The cuvette washer according to Swiegot consists of the following components:

- ① Washing head with function stop (FA30), with lateral screw neck GL14 with cap and suitable hose connection SVS 8 mm **No. 0 185 001, 2190, 111A (2190, 211A)**
- ② Silicone sealing 18 mm, with bore, **No. 2 190 111 B**
- ③ Screw cap GL45 with bore, **No. 0 176 045**
- ④ Silicone/PTFE-sealing 42 x 26 mm, **No. 0 176 026**
- ⑤ Laboratory bottle, GL45, 500 ml, witeg, borosilicate glass 3.3, without cap and pouring ring, **No. 5 527 500 B**
- ⑥ Pouring ring for GL45 bottles (PP), **No. 0 175 202**

Optional:

- 1 vacuum pump Rocker 300, max. vacuum -680 mmHg, 23 l/min, oil and maintenance free, with CE certificate, 230 V, 50/60 Hz, **No. 00829 167300-22**
- 1 m silicone tubing 8 x 2 mm, **No. 9 316 012**
- 1 kg Witonex cleaning concentrate, **No. 5 200 001**

2 990 111

The cuvette washer is usable for cleaning and drying of cuvettes, reaction vessels, small measuring devices like measuring flasks, trapezoidal flasks up to a diameter of approx. 20mm. For safe operation, please wear protective goggles and laboratory gloves.

The components 1-6 must be assembled as follows:

- a. Take the washing head ① and hold it slightly tilted. Insert the screw neck into the bore of the screw cap ③ up to the function stop. Afterwards, put the silicone sealing ② on. Pay attention that the PTFE side is at the bottom meaning that it is directed to the bottle neck.
- b. Screw this component assembly onto the laboratory bottle ⑤, after you put the pouring ring ⑥ in onto the bottle neck.
- c. Laterally, the SVS-screw cap/hose connection unit has to be added to the head ① so that you can later on connect a vacuum pump or respectively a water jet pump.

We recommend our vacuum pump Rocker 300, which we offer optionally including a suitable silicone tubing 8 x 2 mm.

d. Fill now the washing solution into the chamber of the head. As washing solution you can either use your standard liquid cleanser or our laboratory glass cleaning concentrate Witonex Nr. 5 200 001.

The washer can also be used with hot washing solutions (up to +80 °C) without any problem (approx. 100 ml washing solution with 2 990 111 and 2 x 50 ml with 2 990 211)

e. Insert the silicone sealing ② centrally in the spray nozzle channel. The spray nozzle must be well visible.

f. Before you activate the vacuum or water jet pump, tightly press the cuvette (or other glass parts) with the open side onto the silicone tubing.

g. Thereafter, you can activate the pump and repeat the cleaning process (max. 5x) until the dirty article is completely clean and the laboratory bottle is full.

h. The washing solution is collected in the laboratory bottle 500 ml and can possibly at minor impurity be filtered and reused in short time (ex. our respective vacuum filtration system 00829 189302-22).

i. Since you carefully fix each article individually, our cuvette washer offers a gentle purification of glass parts, which can be damaged during automatic cleaning processes.

j. The drying of the respective glass parts is executed similarly without liquid (washing solution).

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