

WITO pipetting tools

Pipette filler with 3 valves

made of natural rubber, standard model for pipettes up to 100 ml, resistant to age and chemicals

	Color	Order number
Standard up to 10 ml	red	3 965 000
Standard up to 10 ml	black	3 965 001
Standard up to 10 ml	green	3 965 002
Standard up to 10 ml	blue	3 965 003
Universal with elongated cone up to 100ml	red	3 966 000
Universal with adapter up to 100ml	red	3 965 100


NEW COLOURS



Flip-Pipette filler with 2 valves

easy to use, two pressure valves for suction and release of liquid up to 100 ml, one automatic air purge valve, effective cleaning by removing the automatic valve, in 5 different colors available


Color	Order number
red	3 967 000
black	3 967 001
yellow	3 967 002
blue	3 967 003
green	3 967 004



Suction bulb

made of PVC, for pipettes with diameter of 6 mm

Size	Suction volume	Ball diameter	Order number
B oval	5 ml	25 mm	3 597 001
0	8 ml	30 mm	3 597 002
2	25 ml	42 mm	3 597 003
3	40 ml	48 mm	3 597 004
5	89 ml	62 mm	3 597 005
7	143 ml	72 mm	3 597 006
9	224 ml	84 mm	3 597 007
11	347 ml	94 mm	3 597 008



Suction bulbs for pipettes 1 - 15 ml

made of silicone

For capacity	Color	Pack size	Order number
approx. 2 ml	transparent	10 pcs.	3 968 000 S
1 ml	neutral	10 pcs.	3 968 001 S
3 ml	neutral	10 pcs.	3 968 002 S
5 ml	neutral	10 pcs.	3 968 005 S
10 ml	neutral	10 pcs.	3 968 010 S



Universal safety pipette bulb

red rubber and silicone pipette bulb for pipettes up to 25 ml
made of synthetic rubber and silicone

NEW

For capacity	Order number
25 ml	3 595 000



Pipette pump with fast release „pi-pump“

made of PP, suitable for pipettes up to 25 ml, suction and dispensing achieved by rotating the thumbwheel
with exception of the 0.2 ml pi-pump, an additional lever releases vacuum to let liquid flow quickly

For capacity	Color-code	Order number
0.2 ml	yellow	7 032 001
2 ml	blue	7 032 002
10 ml	green	7 032 003
25 ml	red	7 032 004
Stand for „pi-pump“		7 032 111

