

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 053 B

**Burettes, 5 ml, straight stopcock, AS, conf. cert.**




**Batch:** 04-22

**Serial:**

**Measuring Nr.** 53930

**Batch size:** 67

| Volume:                | 5,0000 ml   | 3,0000 ml   | 1,0000 ml   |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| <b>Average value:</b>  | 4,9946 ml   | 2,9937 ml   | 1,0008 ml   |
| <b>Deviation:</b>      | 0,0019 ml   | 0,0016 ml   | 0,0020 ml   |
| <b>Precision (CV):</b> | 0,0371 %    | 0,0525 %    | 0,2034 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 7

**Date:** 31.03.2022

**Balance:** Sartorius -  
LE244S-OCE

DE-M