



Certificate of performance

Article number	3 010 053	Burettes 5 ml class A straight. stopc. PTFE plug	
Batch	20-11	Volume	5,0000ml 2,5000ml
Serial number		Tolerance	0,0100ml 0,0100ml
Measuringnr./Batch size	10012 / 6	Measuring volume	5,0035 ml 2,4997 ml
Balance	Sartorius - LE244S-0CE	Deviation	0,0707 % -0,0109 %
Operator	9	Number of measurings	3 3
Date	09.12.2016		



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
 97877 Wertheim
 Germany
 +49 (0) 9342 / 9301-0
 | info@witeg.de

