


# Certificate of performance

Article number	3 010 053	<b>Burettes 5 ml class A straight. stopc. PTFE plug</b>			
Batch	20-11	Volume	5,0000ml	3,0000ml	1,0000ml
Serial number	-	Tolerance	0,0100ml	0,0100ml	0,0100ml
Measuringnr./Batch size	10078 / 1	Average volume	5,0022 ml	3,0026 ml	1,0007 ml
Balance	Sartorius - LE244S-0CE	Deviation	0,0000 ml	0,0000 ml	0,0000 ml
Operator	9	Number of measurings	1	1	1
Date	03.02.2017				



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
 97877 Wertheim  
 Germany  
 +49 (0) 9342 / 9301-0  
 | info@witeg.de

