

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 103 B

**Burettes, 10 ml, str. stopc., amber grad., AS, conf. cert.**




**Batch:** 13-18

**Serial:**

**Measuring Nr.** 20014

**Batch size:** 200

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9966 ml	5,9964 ml	1,9990 ml
Deviation:	0,0060 ml	0,0067 ml	0,0021 ml
Precision (CV):	0,0604 %	0,1110 %	0,1071 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	19.04.2018
Balance:	Sartorius - LE244S-OCE

