




# Certificate of performance

Article number 3 010 103  
 Batch 36-16  
 Serial number  
 Measuringnr./Batch size 10102 / 30  
 Balance Sartorius - LE244S-0CE  
 Operator 9  
 Date 10.02.2017

## Burettes 10 ml class AS clear glass st.stopcock st.-PTFE-plug

Volume	10,0000ml	6,0000ml	2,0000ml
Tolerance	0,0200ml	0,0200ml	0,0200ml
Average volume	9,9965 ml	5,9958 ml	1,9999 ml
Deviation	0,0039 ml	0,0026 ml	0,0024 ml
Number of measurings	5	5	5



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
 97877 Wertheim  
 Germany  
 +49 (0) 9342 / 9301-0  
 | info@witeg.de

