

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 253 B

Burettes, 25 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 47-15

Serial:

Measuring Nr. 46764

Batch size: 192

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0074 ml	15,0081 ml	5,0113 ml
Deviation:	0,0070 ml	0,0080 ml	0,0052 ml
Precision (CV):	0,0278 %	0,0536 %	0,1046 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	17.05.2021
Balance:	Sartorius - LE244S-OCE

