

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 253

Burettes, 25 ml, str. stopcock, PTFE, AS, conf. cert.




Batch: 15-18

Serial:

Measuring Nr.: 24264

Batch size: 46

| Volume: | 25,0000 ml | 15,0000 ml | 5,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0300 ml | ± 0,0300 ml | ± 0,0300 ml |
| Average value: | 24,9922 ml | 14,9742 ml | 4,9829 ml |
| Deviation: | 0,0046 ml | 0,0013 ml | 0,0031 ml |
| Precision (CV): | 0,0184 % | 0,0089 % | 0,0621 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 26.09.2018 |
| Balance: | Sartorius - LE244S-OCE |

DE-M