

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 255

Burettes, 25 ml, lat. needle valve, AS, conf. cert.




Batch: 06-21

Serial:

Measuring Nr. 64650

Batch size: 20

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0053 ml	14,9947 ml	5,0003 ml
Deviation:	0,0026 ml	0,0072 ml	0,0043 ml
Precision (CV):	0,0102 %	0,0482 %	0,0851 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	24.10.2023
Balance:	Sartorius - MSE224S-0CE -DU

