


Certificate of performance

Article number	3 010 502	Burettes 50 ml class AS st.-stopcock solid plug			
Batch	03-17	Volume	50,0000ml	30,0000ml	10,0000ml
Serial number		Tolerance	0,0500ml	0,0500ml	0,0500ml
Measuringnr./Batch size	10083 / 20	Average volume	50,0128 ml	29,9780 ml	9,9920 ml
Balance	Sartorius - LE244S-0CE	Deviation	0,0078 ml	0,0089 ml	0,0060 ml
Operator	9	Number of measurings	4	4	4
Date	09.02.2017				



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH
 Am Bildacker 16
 97877 Wertheim
 Germany
 +49 (0) 9342 / 9301-0
 | info@witeg.de

