

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 503

**Burettes, 50 ml, str. stopcock, PTFE, AS, conf. cert.**




**Batch:** 38-17

**Serial:**

**Measuring Nr.** 33269

**Batch size:** 100

| Volume:                | 50,000 ml   | 30,000 ml   | 10,000 ml   |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| <b>Average value:</b>  | 50,0051 ml  | 29,9766 ml  | 9,9920 ml   |
| <b>Deviation:</b>      | 0,0170 ml   | 0,0098 ml   | 0,0051 ml   |
| <b>Precision (CV):</b> | 0,0341 %    | 0,0326 %    | 0,0512 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 9                         |
| <b>Date:</b>     | 26.09.2019                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

