

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 504

Burettes 50 ml, lat. stopc., glass-plug, AS, conf. cert.




Batch: 38-17

Serial:

Measuring Nr. 15803

Batch size: 4

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0343 ml	29,9881 ml	9,9922 ml
Deviation:	0,0033 ml	0,0019 ml	0,0038 ml
Precision (CV):	0,0065 %	0,0062 %	0,0378 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	21.11.2017
Balance:	Sartorius - LE244S-OCE

DE-M