

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 902

Burettes, 100 ml, str. stopc., solidplug, AS, conf.cert.




Batch: 25-20

Serial:

Measuring Nr. 66419

Batch size: 2

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	99,9667 ml	59,9334 ml	19,9660 ml
Deviation:	0,0084 ml	0,0023 ml	0,0439 ml
Precision (CV):	0,0084 %	0,0039 %	0,2193 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	18.03.2024
Balance:	Sartorius - LE244S-OCE

