

# Certificate of performance



**Article:** 3 010 903

**Burettes, 100 ml, straight stopc., PTFE, AS, conf. cert.**




**Batch:** 38-16

**Serial:**

**Measuring Nr.:** 33780

**Batch size:** 1

|                        |                   |                  |                  |
|------------------------|-------------------|------------------|------------------|
| <b>Volume:</b>         | <b>100,000 ml</b> | <b>60,000 ml</b> | <b>20,000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,1000$ ml   | $\pm 0,1000$ ml  | $\pm 0,1000$ ml  |
| <b>Measured value:</b> | 99,9432 ml        | 59,9470 ml       | 20,0000 ml       |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 7                         |
| <b>Date:</b>     | 25.10.2019                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

