

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 903

Burettes, 100 ml, straight stopc., PTFE, AS, conf. cert.




Batch: 38-16

Serial:

Measuring Nr. 36531

Batch size: 4

| Volume: | 100,000 ml | 60,000 ml | 20,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,1000 ml | ± 0,1000 ml | ± 0,1000 ml |
| Average value: | 99,9632 ml | 59,9273 ml | 19,9873 ml |
| Deviation: | 0,0086 ml | 0,0109 ml | 0,0035 ml |
| Precision (CV): | 0,0086 % | 0,0182 % | 0,0177 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|-----------------------------------|
| Operator: | 7 |
| Date: | 09.03.2020 |
| Balance: | Sartorius - MSE224S-0CE -DU |

