

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 903

**Burettes, 100 ml, straight stopc., PTFE, AS, conf. cert.**




**Batch:** 25-20

**Serial:**

**Measuring Nr.** 64495

**Batch size:** 1

<b>Volume:</b>	<b>100,000 ml</b>	<b>60,000 ml</b>	<b>20,000 ml</b>
<b>Tolerance:</b>	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
<b>Measured value:</b>	100,0579 ml	60,0337 ml	20,0236 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	26.09.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

