




# Certificate of performance

Article number 3 012 053  
 Batch 24-16  
 Serial number -  
 Measuringnr./Batch size 10077 / 2  
 Balance Sartorius - LE244S-0CE  
 Operator 9  
 Date 02.02.2017

## Burettes 50 ml class AS ST-stopcock Schellbach KB

Volume	5,0000ml	3,0000ml	1,0000ml
Tolerance	0,0100ml	0,0100ml	0,0100ml
Average volume	5,0001 ml	2,9958 ml	0,9994 ml
Deviation	0,0001 ml	0,0010 ml	0,0023 ml
Number of measurings	2	2	2



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
 97877 Wertheim  
 Germany  
 +49 (0) 9342 / 9301-0  
 | info@witeg.de

