




Certificate of performance

Article number 3 012 053
 Batch 24-16
 Serial number
 Measuringnr./Batch size 10088 / 3
 Balance Sartorius - LE244S-0CE
 Operator 9
 Date 09.02.2017

Burettes 50 ml class AS ST-stopcock Schellbach KB

Volume	5,0000ml	3,0000ml	1,0000ml
Tolerance	0,0100ml	0,0100ml	0,0100ml
Average volume	4,9973 ml	2,9975 ml	1,0006 ml
Deviation	0,0035 ml	0,0024 ml	0,0019 ml
Number of measurings	3	3	3



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
 97877 Wertheim
 Germany
 +49 (0) 9342 / 9301-0
 | info@witeg.de



ISO 9001
 CERTIFIED

