

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 053

**Burettes, 5 ml, str. stopc., PTFE- plug, AS, conf. cert.**




**Batch:** 19-22

**Serial:**

**Measuring Nr.:** 57484

**Batch size:** 40

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0012 ml	2,9985 ml	1,0004 ml
Deviation:	0,0038 ml	0,0033 ml	0,0011 ml
Precision (CV):	0,0752 %	0,1098 %	0,1111 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

18.07.2022

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

