

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 104

**Burettes, 10 ml, lat. stopcock, AS, conf. cert.**




**Batch:** 36-16

**Serial:** -

**Measuring Nr.:** 14863

**Batch size:** 4

| Volume:         | 10,000 ml   | 6,000 ml    | 2,000 ml    |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value:  | 10,0102 ml  | 6,0020 ml   | 2,0002 ml   |
| Deviation:      | 0,0011 ml   | 0,0003 ml   | 0,0021 ml   |
| Precision (CV): | 0,0105 %    | 0,0058 %    | 0,1056 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

11.10.2017

**Balance:**

Sartorius -  
LE244S-OCE

DE-M