

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 104

**Burettes, 10 ml, lateral stopcock, AS, conf. cert.**




**Batch:** 43-19

**Serial:**

**Measuring Nr.** 37841

**Batch size:** 5

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0108 ml	6,0034 ml	1,9986 ml
Deviation:	0,0037 ml	0,0028 ml	0,0006 ml
Precision (CV):	0,0372 %	0,0462 %	0,0301 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	13.05.2020
Balance:	Sartorius - MSE224S-0CE -DU

