

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 105

Burette, 10 ml, Schellb., lat. Valve, AS, conf. cert.




Batch: 23-21

Serial:

Measuring Nr. 56713

Batch size: 2

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0116 ml	6,0006 ml	1,9997 ml
Deviation:	0,0027 ml	0,0043 ml	0,0030 ml
Precision (CV):	0,0270 %	0,0720 %	0,1490 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	22.06.2022
Balance:	Sartorius - LE244S-OCE

DE-M