

Certificate of performance



Article: 3 012 253 B

Burettes, 25 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 47-15

Serial:

Measuring Nr. 16441

Batch size: 5

| Volume: | 25,000 ml | 15,000 ml | 5,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0300 ml | ± 0,0300 ml | ± 0,0300 ml |
| Average value: | 25,0206 ml | 15,0107 ml | 5,0105 ml |
| Deviation: | 0,0036 ml | 0,0038 ml | 0,0027 ml |
| Precision (CV): | 0,0144 % | 0,0253 % | 0,0545 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 13.12.2017 |
| Balance: | Sartorius - LE244S-OCE |

