




Certificate of performance

Article number 3 012 253
 Batch 03-17
 Serial number
 Measuringnr./Batch size 10128 / 120
 Balance Sartorius - LE244S-0CE
 Operator 9
 Date 23.02.2017

Burettes 25 ml Class AS KB straight stopc. PTFE plug

Volume	25,0000ml	15,0000ml	5,0000ml
Tolerance	0,0300ml	0,0300ml	0,0300ml
Average volume	24,9954 ml	14,9894 ml	4,9993 ml
Deviation	0,0059 ml	0,0075 ml	0,0034 ml
Number of measurings	10	10	10



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
 97877 Wertheim
 Germany
 +49 (0) 9342 / 9301-0
 | info@witeg.de

