

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 253

Burettes, 25 ml, str. stopc., PTFE plug, AS, conf. cert.




Batch: 41-17

Serial:

Measuring Nr.: 20313

Batch size: 10

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9955 ml	14,9984 ml	5,0031 ml
Deviation:	0,0052 ml	0,0119 ml	0,0006 ml
Precision (CV):	0,0208 %	0,0795 %	0,0119 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	03.05.2018
Balance:	Sartorius - LE244S-OCE

