

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 254

Burette, 25 ml, lat. stopc., glass-plug, AS, conf. cert.




Batch: 09-22

Serial:

Measuring Nr. 63873

Batch size: 15

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9824 ml	15,0078 ml	5,0112 ml
Deviation:	0,0034 ml	0,0088 ml	0,0042 ml
Precision (CV):	0,0137 %	0,0584 %	0,0847 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

27.07.2023

Balance:

Sartorius -
LE244S-OCE

DE-M