

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 255

Burettes, 25 ml, lat. needle valve, AS, conf. cert.




Batch: 18-24

Serial:

Measuring Nr. 71240

Batch size: 30

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0093 ml	15,0013 ml	5,0048 ml
Deviation:	0,0078 ml	0,0046 ml	0,0053 ml
Precision (CV):	0,0312 %	0,0304 %	0,1057 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 77

Date: 18.12.2024

Balance: Sartorius -
LE244S-OCE

DE-M