

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 502

**Burettes, 50 ml, str. stopc., glass-plug, AS, conf. cert**




**Batch:** 16-18

**Serial:**

**Measuring Nr.** 23934

**Batch size:** 27

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	49,9952 ml	29,9853 ml	9,9842 ml
<b>Deviation:</b>	0,0078 ml	0,0036 ml	0,0088 ml
<b>Precision (CV):</b>	0,0156 %	0,0121 %	0,0877 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	11.09.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

DE-M