

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 503

Burettes, 50 ml, str. stopc., PTFE plug, AS, conf. cert.




Batch: 30-17

Serial:

Measuring Nr. 14459

Batch size: 247

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0007 ml	30,0077 ml	10,0009 ml
Deviation:	0,0142 ml	0,0141 ml	0,0129 ml
Precision (CV):	0,0284 %	0,0472 %	0,1294 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	04.10.2017
Balance:	Sartorius - LE244S-OCE

