

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 503

Burettes, 50 ml, str. stopc., PTFE plug, AS, conf. cert.




Batch: 25-20

Serial:

Measuring Nr. 39129

Batch size: 200

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0126 ml	30,0067 ml	10,0139 ml
Deviation:	0,0130 ml	0,0186 ml	0,0113 ml
Precision (CV):	0,0261 %	0,0619 %	0,1127 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	21.07.2020
Balance:	Sartorius - LE244S-OCE

