

Certificate of performance



Article: 3 012 504

Burettes, 50 ml, lat. stopc., Schellbach, AS, conf. cert.




Batch: 43-20

Serial:

Measuring Nr.: 42402

Batch size: 30

| Volume: | 50,0000 ml | 30,0000 ml | 10,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 49,9853 ml | 29,9982 ml | 9,9988 ml |
| Deviation: | 0,0172 ml | 0,0189 ml | 0,0084 ml |
| Precision (CV): | 0,0344 % | 0,0629 % | 0,0841 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|-----------------------------------|
| Operator: | 4 |
| Date: | 11.11.2020 |
| Balance: | Sartorius - MSE224S-0CE -DU |

