

Certificate of performance



Article: 3 012 505

Burettes, 50 ml, lat. needle valve., Schellb., AS, conf. c




Batch: 16-18

Serial:

Measuring Nr.: 22371

Batch size: 20

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0004 ml	30,0065 ml	10,0023 ml
Deviation:	0,0135 ml	0,0147 ml	0,0187 ml
Precision (CV):	0,0271 %	0,0489 %	0,1865 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	19.07.2018
Balance:	Sartorius - LE244S-OCE

