

# Certificate of performance



**Article:** 3 012 902

**Burette, 100 ml, str. stopcock, Schellb., AS, conf. cert.**

**Batch:** 43-19




**Serial:**

**Measuring Nr.:** 54710

**Batch size:**

7

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	99,9366 ml	59,9661 ml	20,0052 ml
Deviation:	0,0104 ml	0,0260 ml	0,0046 ml
Precision (CV):	0,0104 %	0,0433 %	0,0228 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

27.04.2022

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU

