

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 903

**Burettes, 100 ml, straight stopc., PTFE, AS, conf. cert.**




**Batch:** 22-22

**Serial:**

**Measuring Nr.** 66449

**Batch size:** 19

Volume:	100,000 ml	60,000 ml	20,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	100,0086 ml	60,0395 ml	19,9936 ml
<b>Deviation:</b>	0,0146 ml	0,0093 ml	0,0227 ml
<b>Precision (CV):</b>	0,0146 %	0,0155 %	0,1135 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	20.03.2024
<b>Balance:</b>	Sartorius - LE244S-OCE

