

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 905

**Burette, 100 ml, lat. stopc., Schellbach, AS, conf. cert.**




**Batch:** 27-20

**Serial:**

**Measuring Nr.** 39089

**Batch size:** 10

<b>Volume:</b>	<b>100,000 ml</b>	<b>60,000 ml</b>	<b>20,000 ml</b>
<b>Tolerance:</b>	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
<b>Average value:</b>	99,9557 ml	60,0039 ml	20,0214 ml
<b>Deviation:</b>	0,0255 ml	0,0282 ml	0,0170 ml
<b>Precision (CV):</b>	0,0255 %	0,0470 %	0,0852 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	16.07.2020
<b>Balance:</b>	Sartorius - LE244S-OCE

