

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 014 053

**Burettes, 5 ml, str. stopc., amber glass, AS, conf. cert.**

**Batch:** 36-16




**Serial:**

**Measuring Nr.** 16008

**Batch size:**

1

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Measured value:	4,9966 ml	2,9938 ml	0,9970 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

27.11.2017

**Balance:**

Sartorius -  
LE244S-OCE

DE-M