

Certificate of performance



Article: 3 014 053

Burettes, 5 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 36-16

Serial:

Measuring Nr.: 31807

Batch size: 1

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Measured value:	4,9940 ml	2,9934 ml	0,9952 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.07.2019
Balance:	Sartorius - LE244S-OCE

