

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 253

Burettes, 25 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 47-19

Serial:

Measuring Nr. 36213

Batch size: 100

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0016 ml	15,0101 ml	4,9928 ml
Deviation:	0,0083 ml	0,0085 ml	0,0074 ml
Precision (CV):	0,0331 %	0,0564 %	0,1483 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	18.02.2020
Balance:	Sartorius - MSE224S-0CE -DU

