

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 503

Burettes, 50 ml, str. stopc., amber glass, AS, conf. cert.




Batch: 38-17

Serial:

Measuring Nr. 34419

Batch size: 99

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0118 ml	30,0002 ml	10,0064 ml
Deviation:	0,0171 ml	0,0204 ml	0,0146 ml
Precision (CV):	0,0342 %	0,0680 %	0,1463 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.11.2019
Balance:	Sartorius - MSE224S-0CE -DU

