

# Certificate of performance

**Article:** 3 014 505

**Burettes, 50 ml, lat. stopc., amber glass, AS, conf. cert.**




**Batch:** 03-17

**Serial:**

**Measuring Nr.:** 38830

**Batch size:** 9

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0159 ml	29,9821 ml	9,9832 ml
Deviation:	0,0126 ml	0,0112 ml	0,0079 ml
Precision (CV):	0,0253 %	0,0372 %	0,0791 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

25.06.2020

**Balance:**

Sartorius -  
LE244S-OCE