

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 020 257

Burettes, 25 ml, str. needle valve, AS, conf. cert.




Batch: 06-21

Serial:

Measuring Nr.: 60501

Batch size: 20

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0104 ml	15,0002 ml	5,0053 ml
Deviation:	0,0038 ml	0,0009 ml	0,0039 ml
Precision (CV):	0,0154 %	0,0060 %	0,0781 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	28.11.2022
Balance:	Sartorius - MSE224S-0CE -DU

