

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 020 507 DV

**Burettes, 50 ml, 2 needle valve, AS, conf. cert.**




**Batch:** 03-17

**Serial:**

**Measuring Nr.** 14228

**Batch size:** 10

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0221 ml	29,9788 ml	10,0048 ml
Deviation:	0,0062 ml	0,0058 ml	0,0041 ml
Precision (CV):	0,0124 %	0,0194 %	0,0414 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	25.09.2017
Balance:	Sartorius - LE244S-OCE

