

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 020 507

Burettes, 50 ml, straight needle valve, AS, conf. cert.




Batch: 38-17

Serial:

Measuring Nr. 35122

Batch size: 30

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0218 ml	29,9922 ml	9,9914 ml
Deviation:	0,0131 ml	0,0143 ml	0,0035 ml
Precision (CV):	0,0261 %	0,0475 %	0,0353 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.12.2019
Balance:	Sartorius - MSE224S-0CE -DU

