

Certificate of performance



Article: 3 020 907

Burettes, 100 ml, str. needle valve, AS, conf. cert.




Batch: 38-16

Serial:

Measuring Nr. 16841

Batch size: 5

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
Average value:	99,9617 ml	59,9577 ml	19,9836 ml
Deviation:	0,0299 ml	0,0274 ml	0,0367 ml
Precision (CV):	0,0299 %	0,0456 %	0,1834 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	08.01.2018
Balance:	Sartorius - LE244S-OCE

