

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 022 507

Burettes, 50 ml, straight needle valve, AS, conf. cert.




Batch: 01-19

Serial:

Measuring Nr.: 40004

Batch size: 17

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0098 ml	30,0226 ml	10,0273 ml
Deviation:	0,0184 ml	0,0047 ml	0,0062 ml
Precision (CV):	0,0368 %	0,0156 %	0,0621 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	09.09.2020
Balance:	Sartorius - LE244S-OCE

DE-M