

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 022 907

**Burettes, 100 ml, straight needle valve, AS, conf. cert.**




**Batch:** 26-20

**Serial:**

**Measuring Nr.:** 39196

**Batch size:** 15

Volume:	100,000 ml	60,000 ml	20,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	99,9946 ml	59,9962 ml	20,0202 ml
Deviation:	0,0259 ml	0,0101 ml	0,0134 ml
Precision (CV):	0,0259 %	0,0168 %	0,0672 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	24.07.2020
Balance:	Sartorius - LE244S-OCE

