

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 024 257

**Burette, 25 ml, str. needle Valve, amber, AS, conf. cert.**




**Batch:** 47-19

**Serial:**

**Measuring Nr.** 38009

**Batch size:** 2

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0190 ml	14,9867 ml	5,0023 ml
Deviation:	0,0038 ml	0,0058 ml	0,0166 ml
Precision (CV):	0,0152 %	0,0388 %	0,3312 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	15.05.2020
Balance:	Sartorius - LE244S-OCE

DE-M