

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 024 507

Burette, 50 ml, str. needle valve, amber, AS, conf. cert.




Batch: 03-17

Serial:

Measuring Nr. 34653

Batch size: 2

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9998 ml	29,9723 ml	10,0043 ml
Deviation:	0,0013 ml	0,0074 ml	0,0098 ml
Precision (CV):	0,0027 %	0,0246 %	0,0984 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	26.11.2019
Balance:	Sartorius - MSE224S-0CE -DU

DE-M