

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 105

**Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 43-19

**Serial:**

**Measuring Nr.** 38467

**Batch size:** 10

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0024 ml	6,0000 ml	1,9973 ml
Deviation:	0,0061 ml	0,0038 ml	0,0078 ml
Precision (CV):	0,0605 %	0,0638 %	0,3913 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	09.06.2020
Balance:	Sartorius - LE244S-OCE

